



## **RED**30-S

No other stereo microscope stand can be compared either in design or quality. The base is completely sealed to protect the illuminator from liquids that might be spilled. All the models are fully equipped with incident and transmitted LED illumination, carrying handle for easy transportation and "one touch" sample clips.

Model RED 30-S

Optical System Greenough

Observation Tube Binocular head

Interpupillary distance 52-69mm

**Diopter adjustment** On the left tube, +/- 5 diopter

Inclination 45° inclined
Eyepieces WF10X/20mm

**Objectives system** Turret changer with LED incident illumination

Objectives magnification 2X - 4X

Working distance 50mm

**Stand** Fixed arm type stand

Base 225x179mm Column/Arm 200mm high

Head holder Fixed

Focus mechanism Coarse focusing system with tension adjustment

Focusing stoke 47mm

Incident illuminationLED 3.5V/100mW with intensity controlTransmitted illuminationLED 3.5V/100mW with intensity control

Transformer External

**Power supply** Charger 100-240V or rechargeable batteries

**Accessories included** Black, white and frosted stage plates, dust cover, charger and batteries

Length x Width x Height 225x179x250mm

Weight 3.5Kg



Add a Moticam camera to capture, document, annotate and share images and videos with the Motic Images Plus 2.0 software.